


# A Beginners Guide to Phred, Phrap and Consed

Training offered by the Advanced  
Biomedical Computing Center,  
NCI/Frederick


Beena Neelam  
email: [bneelam@ncifcrf.gov](mailto:bneelam@ncifcrf.gov)  
call: 301 846 5779



# Objectives


- **Objectives of the course**

- What are phred, phrap and consed?
- prepare your chromat files for running phred and phrap on the ABCC UNIX machines
- run phred- base calling program
- run phrap- assembly program
- view contigs using consed- a visualization tool



# What are Phred/Phrap/Consed?

- **Phred/Phrap/Consed is a distributed package, free for academic use, for:**
  - Reading trace files (chromatograms)
  - Quality value assignment to each base
  - Vector and repeat sequences masking
  - Sequence assembly
  - Assembly visualization and editing
  - Automatic finishing



# Phred


- **Developed at the University of Washington**

- Brent Ewing, et. al., Base-calling of automated sequencer traces using phred. I. Accuracy assessment. 1998. Genome Research 8:175-185.
  - Brent Ewing and Phil Green, Base-calling of automated sequencer traces using phred. II. Error probabilities. 1998. Genome Research 8:186-194.
- 
- Reads DNA sequencer trace data from chromatogram in SCF (standard chromatogram format), ABI (373/377/3700), ESD (MegaBACE) and LI-COR
  - calls bases (sequencing machine and chemistry specific) contd.




# Trace File Quality


High



Medium




Low



# Phred

- assigns quality score to bases (related to base call error probability)
- Determines file format automatically
- Writes sequences to output in FASTA format or PHD format or SCF (Sequence Comparison Interchange Format) format
- Outputs quality scores for each trace file
- Optionally outputs Polyphred values file




# phd2fasta

- **phd2fasta program**
  - converts .phd files to sequence in multifasta format
  - writes .qual file (quality scores) for each trace file
- **Output:**
  - fasta.seq contains fasta sequences
  - fasta.seq.qual contains quality scores



# Crossmatch

- A general purpose utility for comparing two DNA sequences-efficient Smith Waterman algorithm
- More sensitive than BLASTN but slower, allows gaps.
- Can be used to
  - compare a set of DNA sequences to a set of vector sequences
  - compare contigs derived by two methods of assembly
  - align contigs to reference sequence



# Phrap

- **F(ph)ragment Assembly Program- Phil Green**
  - A program for assembling DNA sequence data-reads of 500-700 bp,
  - Uses entire read, uses user supplied and internally computed quality values
  - Constructs contigs as a mosaic of highest quality reads not as consensus



# Phrap

- gives extensive information about the assembly -easier to troubleshoot
  - in <filename>.phrap.out file
    - Residue counts, Read name analysis, Input quality, Regions converted to N's
- can handle large datasets



# Consed

- **A program for viewing and editing sequence assemblies -David Gordon**

- Gordon, D., C. Desmarais, and P. Green. 2001. Automated Finishing with Autofinish. Genome Research. 11(4):614-625.
- The version we are using is consed13.0
- View contigs assembled with Phrap, view single or multiple trace
- Add reads and reassemble
- Several useful and timesaving navigation features to view aligned reads, low/high quality regions
- Autofinish: Close gaps, improve sequence quality, determine the relative orientation of contigs
- several other features...



# Running the programs

- **Two ways to run the programs**
  - Run as a single program called PhredPhrap
    - PhredPhrap program runs Phred, runs phd2fasta to make the fasta format sequence file from the .phd files, runs cross match to mask vector and then runs phrap to assemble the reads
  - Run sequentially
    - Run Phred, create .phd files, then run phd2fasta (providing quality scores if necessary), then run crossmatch to mask vector and then run phrap



# Running Phred (base calling)

- **Make directories**


- chromat\_dir, phd\_dir and edit\_dir and place chromatogram files in chromat\_dir

- **Run phred**

- either compressed file or all the individual trace files
- phd\_dir will have .phd files


- **Run phd2fasta**

- with .phd files as input, and edit\_dir as output dir. Result is a set of fasta files in edit\_dir



# phredpar.dat file


- **phredpar.dat file**
  - is a chemistry parameter file, contains dye primer data
- **Location of phredpar.dat file:**
  - ncisgi:            /abcc/apps/IRIX64/phred02/



- **Example of a phredpar.dat file:**
- #####
- # phredpar.dat - phred parameter file: 020425
- # known chemistries: primer, terminator, unknown
- # known dyes : rhodamine, d-rhodamine, big-dye energy-transfer, bodipy, unknown
- # known machines : ABI\_373\_377, MolDyn\_MegaBACE, ABI\_3700, LI-COR\_4000,
- # Beckman\_CEQ\_2000, ABI\_3100
- # Notes:
- # (1) enclose the `dye primer` name in double quotes and include spaces in the names.
- # (2) leave one or more spaces between the `dye primer` and chemistry names, between
- # the chemistry and dye names, and between the dye and machine names.
- # (3) add entries between the `begin chem\_list` and `end chem\_list` lines.
- # (4) `dye primer` string `\_\_phred\_default\_\_` identifies the chemistry, dye, and machine
- # used when the
- # chromatogram dye primer string is missing from either the chromatogram or this file and
- # phred is run with the -process\_nomatch option.
- #####
- begin chem\_list
- "" primer big-dye ABI\_3700
- "DP6%25Ac{-21M13}" primer rhodamine ABI\_373\_377
- "DP6%Ac{-21M13}" primer rhodamine ABI\_373\_377
- "DP6%25Ac{M13Rev}" primer rhodamine ABI\_373\_377

# Setting up to run the programs

- **ABCC UNIX account**
  - <http://www2.ncifcrf.gov/>
- **Transferring trace files**
  - Finch
  - ftp from desktop
  - as compressed file on CD
  - already on UNIX machine
- **Edit your .cshrc file:**
  - add line `source /abcc/apps/perl/setup.csh`



# phred, phrap and crossmatch location

## Phred on ncisgi:

/abcc/apps/IRIX64/phred02/phred

## Phrap on ncisgi:

/abcc/apps/IRIX64/phrap02/phrap


## cross\_match on ncisgi:

/abcc/apps/IRIX64/phrap02/cross\_match



# Preparing trace files for analysis

- `%cd tar_folder`
- `%gunzip ab1_F.tar.gz`
- `%tar -xf ab1_F.tar`
- `-xf`, extract (x) files (f) from an archive
- Result: list of .ab1 files



# Running phred

- **Make subdirectories**

- in the directory containing the compressed file:
- make chromat\_dir phd\_dir edit\_dir

- **Run phred**

- %phred -id chromat\_dir -pd phd\_dir

- **Output**

- chromat\_dir has the chromat files eg \*.ab1 format
- phd\_dir will have corresponding \*.phd files created by phred



# Phred options

- **Command line options**

- options are preceded by a '-'
- input options, processing options and output options
- input options: -id, -if
- processing options: -exit\_nomatch
- output options: -sa, -qa, -sd, -qd, -pd



# Running phd2fasta

- **phd2fasta program**
  - converts .phd files to sequence in multifasta format
  - writes .qual file (quality scores) for each trace file
  - %phd2fasta -id phd\_dir -os edit\_dir/fasta.seq -oq edit\_dir/fasta.seq.qual
- **Output:**
  - fasta.seq contains fasta sequences
  - fasta.seq.qual contains quality scores



# An example of a .phd file

- BEGIN\_SEQUENCE 522e10\_R\_N1-9H12.g
- BEGIN\_COMMENT
- CHROMAT\_FILE: 522e10\_R\_N1-9H12.g
- ABI\_THUMBPRINT: 0
- PHRED\_VERSION: 0.020425.c
- CALL\_METHOD: phred
- QUALITY\_LEVELS: 99
- TIME: Fri Jan 23 11:41:30 2004
- TRACE\_ARRAY\_MIN\_INDEX: 0
- TRACE\_ARRAY\_MAX\_INDEX: 14306
- TRIM: 23 878 0.0500
- TRACE\_PEAK\_AREA\_RATIO: 0.0091
- CHEM: term
- DYE: big
- END\_COMMENT
- BEGIN\_DNA
- g 6 7
- g 6 16
- a 6 37
- t 6 40
- t 8 60
- g 3 71
- c 2 80



# Running crossmatch to screen vector

- **Run crossmatch**

- %cross fasta.seq  
/abcc/apps/IRIX64/phrap02/allvector.seq  
–minmatch 10 –minscore 20 –screen > screen.out

- **The –screen option**

- causes files called fasta.seq.screen to be created, contains vector masked sequences.

- **screen.out**

- information about screen parameters and vector matches found



# Running phrap

- phrap = F(ph)ragment assembly program or Phil's revised assemble program
- %phrap -new\_ace edit\_dir/fasta.seq > phrap.out
- All the seqs. to be assembled are typically given in a single input file.



- Input files for Phrap: Sequence file (typically single multifasta file) and Quality file, both fasta format
- for each input fasta.seq file if a corresponding fasta file containing data quality information is provided, it greatly improves the accuracy of assembly.
- if C05D11.reads.screen is your input file corresponding quality file should be C05D11.reads.screen.qual




# Phrap options

- Scoring of pair wise alignments
- Banded search
- Filtering of matches
- Input data interpretation
- Assembly
- Consensus sequence construction
- Output

# Phrap output files


- **\*.contigs – fasta file containing the contigs**
  - Contigs with more than one read
  - Singletons (single reads with a match to some other contig but that couldn't be merged consistently to it)
- **\*.singlets – fasta file of the singlet reads**
  - Reads with no match to other read
- **\*.ace**
  - allows for viewing the assembly using Consed
- **\*.view**
  - required for viewing the assembly using Phrapview



# Running Consed

- Viewing and editing sequence assemblies with Consed
  - %cd edit\_dir
  - Type consed13

# Consed windows



**Consed Main Window**

File Navigate Info Options Help

2.ace

Search for String Undo Edit... Print Quit Consed

Assembly View Add New Reads Digests Remove Reads Miniassembly

Single Click and then  
Double Click to Sel

Contig List

Contig1 (17 reads, 15

Read List

2-G3\_A Contig1 (~94  
2-G3\_B Contig1 (~92  
2-G3\_C Contig1 (~88  
2-G3\_D Contig1 (~91  
2-G3\_E Contig1 (~92  
2-G3\_F Contig1 (~92  
2-G4\_C Contig1 (~91

Find reads containing

Find 1st read startin

**Aligned Reads**

File Navigate Info Color Din Misc Help

2.ace Contig1 Sone Tags Pos: I clear

Search for String Compl Cont Compare Cont Find Main Win Err/10kb: 20.72

760 770 780 790 800 810

CONSENSUS GGCCATTGACAGAAGAAAAAATAAAAGCATTAGTAGAAATTTGTACAGAAATGGAAAAGG

hiv-p6pront.ref.scf GGCCATTGACAGAAGAAAAAATAAAAGCATTAGTAGAAATTTGTACAGAAATGGAAAAGG

2-G7\_A GGCCATTGACAGAagaAAAAAATAAAAGCATTAGtagAAATTTGTACAGAAATGGAAAagg

2-G3\_A GGCCATTGAcagaagaAAAAAATAAAAGCATTAGTAGAAATTTGTACAGAAATGGAAAagg

2-G7\_B GGCCATTGACAGAAGAAAAAATAAAAGCATTAGTAGAAATTTGTACAGAAATGGAAAAGG

2-G3\_B GGCCATTGACAGAAGAAAAAATAAAAGCATTAGTAGAAATTTGTACAGAAATGGAAAAGG

2-G7\_E GGCCATTGACAGAAGAAAAAATAAAAGCATTAGTAGAAATTTGTACAGAAATGGAAAAGG

2-G3\_E GGCCATTGACAGAAGAAAAAATAAAAGCATTAGTAGAAATTTGTACAGAAATGGAAAAGG

2-G4\_E GGCCATTGACAGAAGAAAAAATAAAAGCATTAGTAGAAattGTACAGAAATGGAAAAGG

2-G7\_F ggcatttgacagaaGaaaaaATAAAagcattagtagaaAtttgtacagaAATGGAAAagg

2-G3\_F GGCCtttgacagaaGaaAAAAATaaAagcattaATAGAAatttgtacagAAATGGAAAAGg

2-G4\_F gGcctttgacagAAGAAAAAATAAAAGcattaATAGaaaacttgtacAGAAATGGAAAAGg

<<< << < Prev Next > >> >>> cursor dismiss

Show Contig Close All Windows



# X Aligned Reads

File Navigate Info Color Dim Misc

2.ace Contig1

Search for String Compl Cont Compare Cont Find Main Min Err/10

770 780 790

CONSENSUS

hiv-p6prort.ref.scf

2-G7\_A

2-G3\_A

2-G7\_B

2-G3\_B

2-G7\_E

2-G3\_E

2-G4\_E

2-G7\_F

2-G3\_F

2-G4\_F

tag types: dataNeeded tag data

## X Trace Window: Contig1

Dismiss

2-G3\_B con rd

H con

H edt

V phd

Scroll Together? Yes No

Remove Undo

2-G7\_E con rd

H con

H edt

V phd

Scroll Together? Yes No

Remove Undo

2-G7\_B con rd

H con

H edt

V phd

Scroll Together? Yes No

Remove Undo

Help Insert prev Dismiss next Help Delete